

Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.



Iron

Alternative Name

Prussian Blue Iron, Ferric Iron

Methodology

Immunohistochemistry (IHC)

Test Description

Special stain. This iron stain is used to detect and identify ferric (Fe3+) iron in tissue preparations, blood smears, or bone marrow smears. Minute amounts of ferric acid are commonly found in bone marrow and in the spleen. Abnormal amounts of iron can indicate hemochromatosis and hemosiderosis. Lack of iron identification can indicate iron-deficient anemia.

Specimen Requirements

- A formalin-fixed, paraffin-embedded (FFPE) tissue block is preferred specimen type
- or
- One (1) unbaked, unstained slide for H&E staining (required) and two to three (2-3) positively charged unstained slides (all cut at 4-5 microns) for each test/antibody ordered
- Block and slide identifiers should be clearly written and match exactly with the specimen ID and specimen labeling as noted on the requisition.

Storage & Transportation

Use cold pack for transport. Make sure cold pack is not in direct contact with specimen.

CPT Code(s)*

88313x1

New York Approved

Yes

Level of Service

Stain Only

Turnaround Time

24 hours

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way Fort Myers, FL 33912 Phone: 239.768.0600/ Fax: 239.690.4237 neogenomics.com © 2024 NeoGenomics Laboratories, Inc. All Rights Reserved. All other trademarks are the property of their respective owners Rev. 042824